

Study abstract demonstrates utility of Entegriion VCM system for anticoagulation monitoring of COVID-19 patients

Novel VCM showed correlation between clotting time and unfractionated heparin (UFH) superior to “gold standard” aPTT method in single-center study abstract

DURHAM, NC, USA, September 29, 2020 /EINPresswire.com/ -- Entegriion, Inc., a leading developer of innovative solutions for clinical diagnosis and

management of hemostasis, highlighted a newly released abstract of a single-center study demonstrating the clinical value of the company's portable Viscoelastic Coagulation Monitor (VCM) in [detecting hypercoagulation among COVID-19 patients](#). The abstract reported that



“

Providing clinicians with diagnostic tools that simplify the monitoring of critically ill patients is crucial to improving outcomes during the COVID-19 pandemic”

*John B. Mowell, Entegriion
Chairman and Executive
Director*

unfractionated heparin (UFH) was better correlated with the VCM than with Activated Partial Thromboplastin (aPTT), considered the “gold standard” for UFH monitoring. The Entegriion VCM has a CE mark for use in Europe; it is not cleared for use in the U.S.

In the study, conducted at Fondazione IRCCS Ca' Granda Ospedale Policlinico (Milan, Italy) from March 25 to June 8, 2020 under the direction of Principal Investigator Mauro Panigada, VCM was used as a daily assessment of the viscoelastic properties of 36 [critically ill COVID-19 patients](#) and the results compared with standard aPTT coagulation


tests. Anticoagulation was provided with either UFH (4 patients, 52 samples) or with LMWH (low molecular weight heparin) (31 patients, 86 samples). UFH influenced all VCM parameters toward hypocoagulability; LMWH did not. UFH dose was better correlated with VCM clotting time (CT) than with aPTT.

“Providing clinicians with diagnostic tools that simplify the monitoring of critically ill patients is crucial to improving outcomes during the COVID-19 pandemic,” said John B. Mowell, Entegriion Chairman and Executive Director. “The recent study in Milan showed how the Entegriion VCM can

play an important role in helping support clinical decision-making during the COVID-19 pandemic by providing an accurate, easy-to-use device for [monitoring anticoagulation therapy at the patient's bedside.](#)"

About the Viscoelastic Coagulation Monitor (VCM)

The Entegriion VCM is a compact, portable device that performs a viscoelastic analysis of the coagulation-fibrinolysis process utilizing glass surface activation of untreated whole blood. The wide surface area of contact between the blood and the glass inside the cartridge accounts for the rapid initiation of clotting, eliminating the need for activating factors. VCM has been compared to ROTEM® NATEM (non-activated method), showing good-to-moderate agreement in test results between the two systems in a cohort of patients undergoing major surgery.¹ By making viscoelastic measurements of the hemostasis of blood samples rapidly and accurately, the VCM system addresses the accessibility, robustness and training issues associated with large, complex, traditional systems.



VCM Analyzer

About Entegriion

Entegriion is a life sciences company providing innovative solutions for clinical diagnosis and management of hemostasis. Our products bring about major changes to the critical pathways in the diagnosis and treatment of patient coagulopathies, enabling clinicians to rapidly and effectively advance care in both the hospital and field environments.

□For more information, contact info@entegriion-vcm.com

1. Brearton C, Rushton A, Parker J, Martin H, Hodgson J. Performance Evaluation of a New Point of Care Viscoelastic Coagulation Monitoring System in Major Abdominal, Orthopaedic and Vascular Surgery. Platelets. 2020;1–8.

Bert Valada

Entegriion, Inc.

+1 919-536-1500

[email us here](#)

This press release can be viewed online at: <https://www.einpresswire.com/article/527304216>

EIN Presswire's priority is source transparency. We do not allow opaque clients, and our editors try to be careful about weeding out false and misleading content. As a user, if you see something we have missed, please do bring it to our attention. Your help is welcome. EIN Presswire,

Everyone's Internet News Presswire™, tries to define some of the boundaries that are reasonable in today's world. Please see our Editorial Guidelines for more information.

© 1995-2020 IPD Group, Inc. All Right Reserved.